

# Paediatric Elastic Wrist **Thumb Brace**

Material - Elastic



Manufactured  
in the UK



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



Breathable



ISO:13485  
Accredited

Manufactured in the UK; the Elastic Wrist Thumb brace is made from elastic, which allows for gentle compression. The wrist brace gives additional support to the first CMC and MCP joint and has a removable stay can be shaped. The palmar bar and thumb extension bar can be shaped, to provide a more custom fit. The brace is coloured to size with Blue - PS1, Red -PS2 and Green - PS3, for alternative colours contact customer services.



## Indications for use

- » Healing Fractures
- » Healing ligament injuries
- » Juvenile RA
- » Strains and sprains

## Features and benefits

- » Wrist Brace with palmar bar and thumb extension bar which can be shaped to provided a more custom fit
- » Additional support to the first CMC and MPC joints
- » Can be hand washed and air dried
- » Different colours available on request, contact customer services

### Sizing Information

Size	Wrist Circ (mm)	Length (inches)	Colour	Code Right	Code Left
PS1	95 - 115	6	Blue / Black	EWTB/1R/PUP	EWTB/1L/PUP
PS2	115 - 130	6	Red / Black	EWTB/2R/PUP	EWTB/2L/PUP
PS3	130 - 150	6	Green / Black	EWTB/3R/PUP	EWTB/3L/PUP



**Beagle**  
ORTHOPAEDIC