

Paediatric Standard Wrist Brace

Material - Elastic



Manufactured in the UK



Hand Wash & Air Dry Only



Same Day Dispatch



Breathable



ISO:13485 Accredited

Manufactured in the UK and made from elastic, the Paediatric Standard Wrist Brace provides support to the whole wrist region for added reassurance. The wrist brace maintains full finger movement for rehab exercises and has an easy to modify metal stay that keeps the wrist in slight extension, relieving strain on the tendons. The brace is coloured to size with Blue - PS1, Red -PS2 and Green - PS3, for alternative colours contact customer services.



Indications for use

- » Healing Fractures
- » Strains and sprains
- » Juvenile RA

Features and benefits

- » Elastic provides support to the whole wrist region for added reassurance
- » Maintains full finger movement for rehabilitation exercises
- » Easy to modify metal stay keeps wrist in slight extension relieving strain on the tendons
- » Can be hand washed and air dried
- » Different colours available on request, contact customer services

Sizing Information

Size	Wrist Circ (mm)	Length (inches)	Colour	Code Right	Code Left
PS1	95 - 115	6	Blue	SWB/1R/PUP	SWB/1L/PUP
PS2	115 - 130	6	Red	SWB/2R/PUP	SWB/2L/PUP
PS3	130 - 150	6	Green	SWB/3R/PUP	SWB/3L/PUP



Beagle
ORTHOPAEDIC