

# BioSkin OA Spiral Brace


**Material** - Ultima



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



Breathable



ISO:13485  
Accredited

The OA Spiral is a lightweight, low-profile, knee brace for tibiofemoral and/or patellofemoral OA pain. Tibiofemoral: The OA Spiral uses a three-point leverage system and a non-elastic strap to pull the knee out of varus or valgus and decrease pain. Patellofemoral: Among patients with tibiofemoral OA, 91% also have radiographic patellofemoral OA. BioSkin's Q strap has been clinically proven to decrease pain in the PF joint by supporting proper patellar alignment and increasing joint surface contact. The lightweight hinge and thin breathable materials yield a product that is comfortable and easy to wear, even under clothing. The strap and buckle design makes donning and doffing easy for patients who struggle with dexterity. The OA Spiral is also designed for longevity with easily replaceable hook & loop panels so the patient can get years of use from the same brace.



## Indications for use

- » Osteoarthritis
- » Conservative/postoperative
- » Pre-operative
- » Knee arthritis
- » Bone bruising

## Features and benefits

- » Clinically validated pain relief
- » High patient satisfaction
- » Enhanced pain relief
- » Adjustable sizing
- » Contoured cut
- » Improved proprioception
- » Biomechanical improvements
- » Easy to fit, apply and use
- » Hand wash and air dry

# Lateral

## Sizing Information

### Circumference with knee fully extended

Size	Mid - patella	180mm above	Left Code	Right Code
Small	330 - 360	460 - 510	51501L-LAT	51501R-LAT
Medium	360 - 410	510 - 560	51502L-LAT	51502R-LAT
Large	410 - 430	560 - 640	51503L-LAT	51503R-LAT
XLarge	430 - 460	640 - 710	51504L-LAT	51504R-LAT

Contact customer service for other sizes.

# Medial

## Sizing Information

### Circumference with knee fully extended

Size	Mid - patella	180mm above	Left Code	Right Code
Small	330 - 360	460 - 510	51501L-MED	51501R-MED
Medium	360 - 410	510 - 560	51502L-MED	51502R-MED
Large	410 - 430	560 - 640	51503L-MED	51503R-MED
XLarge	430 - 460	640 - 710	51504L-MED	51504R-MED

Contact customer service for other sizes.