

BioSkin Weil Osteotomy Strap


Material - Latex Free



Hand Wash &
Air Dry Only



Same Day
Dispatch



Breathable



ISO:13485
Accredited

The Weil Osteotomy Strap with Compression Wrap, promotes proper healing after a Weil Osteotomy procedure and prevents dorsal drift, floating toe, and hammertoe. The strap wraps comfortably around the toe and anchors it into a neutral position. The entire strap is trimmable for a customized fit and can easily be worn inside a walking boot, post-op shoe, or normal footwear. The pack includes 2 Weil Osteotomy Straps and 1 Compression Wrap.



Indications for use

- » Hammer toe
- » Floating toe
- » Dorsal drift
- » Post Weil Osteotomy surgery

Features and benefits

- » Promotes proper toe alignment
- » Universal sizing
- » Use with AFTR braces or the mid-foot compression wrap
- » Trimmable for a custom fit
- » Comfortable in shoes
- » Perfect for post surgery

Pack includes

- » 2 x Weil Osteotomy Straps
- » 1 x CMC Wrap

Sizing Information

Size	Code Single
Universal	953713