

# Orthomerica California Compression Orthosis



Hand Wash &  
Air Dry Only



Same Day  
Dispatch



ISO:13485  
Accredited

The California compression Orthosis is a low-profile lumbar orthosis that provides maximum abdominal support with one easy pull. Its one-handed compression pull tab makes application simple and quick.

The breathable, comfortable mesh fabric minimises heat buildup common with other spinal orthoses.



## Indications for use

- » Acute and chronic low back pain
- » Spondylolisthesis
- » Post-op low lumbar laminectomy
- » Lightweight construction is optimal for patients
- » Requiring support for everyday activities such as lifting, long periods of standing, work-related motions, golf and gardening

## Features and benefits

- » Low profile lumbar orthosis providing maximum abdominal support with one easy pull
- » Simple and quick
- » Breathable & comfortable mesh fabric minimises heat build up
- » Can be hand washed and air dried



### Sizing Information

Size	Waist Circ (mm)	8" depth Code	10" depth Code
XSmall	635 - 762	B/3199	B/3209
Small	763 - 889	B/3200	B/3210
Medium	890 - 1016	B/3201	B/3211
Large	1017 - 1143	B/3202	B/3212
XLarge	1144 - 1270	B/3203	B/3213

**Beagle**  
ORTHOPAEDIC