

Orthomerica California

ECO Brace

Simple & easy way to stabilise
the pelvis & lumbar spine to ease pain



**Simple to apply as a
backpack, with easy-to-
reach self-adjusting &
padded shoulder straps.**

**The metal
spinal frame
can be easily
contoured to
the patient's
thoracic &
lumbar spine.**

**CONTACT US
TODAY**

E: info@beagleortho.com

T: 01254 268 788

W: www.beagleorthopaedic.com

Beagle[®]
ORTHOPAEDIC

The perfect

SOLUTION

- » Thoracic compression fractures
- » Kyphosis
- » Early through late stages of osteoporosis
- » Strain in upper throacic spine due to sustained positioning

Designed with

FEATURES

- » Contourable, aluminium “spine” is lined with CoolFoam & extends from T2 - S1
- » Patented California™ Compression abdominal component, made of cool, breathable mesh provides cinching mechanical advantage
- » Detachable abdominal component can be used independently of remaining system to provide a functional lumbar spinal orthosis
- » Durable, one-pull compound closure provides patient adjustable compression
- » Patient controlled shoulder extension improves compliance & tolerance
- » Non-twisting strap design makes application easy for elderly individuals with limited shoulder mobility
- » Requires no tools for fitting procedure
- » Padded shoulder straps loosen without completely unfastening in order to facilitate simple reapplication

Sizing Information *Tall usually recommended for those 5'9" & taller

Size	Waist Circ (mm)	Standard Code	Tall* Code
XSmall	635 - 760	B/3253	B/3273
Small	760 - 890	B/3254	B/3274
Medium	890 - 1010	B/3255	B/3275
Large	1010 - 1140	B/3256	B/3276
XLarge	1140 - 1270	B/3257	B/3277