

Orlando® HKAFO (Hip - Knee-Ankle-Foot) Femoral Fx Orthosis Custom-to-Measurement Product

Measurement Legend


- A** - Waist Circumference
- B** - Hip Circumference
- C** - Proximal Thigh Circumference
- D** - Distal Thigh Circumference
- E** - Largest Calf Circumference
- F** - Circumference just above Ankle
- G** - Length, Trochanter to Knee
- H** - Anatomical Length, Perineum to floor
- I** - Caliper Measurement, Knee Centre (ML)
- J** - Length, Knee Centre to floor

Dimensions Unit - Please Tick

Inches Centimetres Millimetres

Left Right

Brace supplied with two femoral socks



Component Options (Circle Choice)

Knee Joint

Shoe Inserts



Free Motion



Drop Lock 1/8" or 1/16"



Adjustable Motion

'Virtual' Hip joint as standard



Free Motion



Solid Ankle

Patient Information

Last Name.....
 Approx' Age.....
 Sex.....
 Height.....
 Weight.....
 Other Info/Instructions.....

Order Information

Date.....
 Order No.....
 Company/Hospital.....
 Address.....
