

Orthomerica Paediatric AFO



Hand Wash &
Air Dry Only



ISO:13485
Accredited

The Orthomerica paediatric AFO keeps a child's feet and ankles in a good position for standing and walking. The AFO is designed to prevent equinus contractures and aid in the treatment of patients with spastic hemiplegia or diplegia.

The low-profile design is prefabricated and made from polypropylene that can easily be fitted into a child's shoes. The leaf spring design allows for a smoother gait pattern and low-cost blanks can be custom fitted. Especially effective application in clinical and office settings.



Indications for use

- » Flaccid foot drop
- » Equinus contractures
- » Spastic hemiplegia
- » Diplegia

Features and benefits

- » Popular prefabricated design accommodates many patients not requiring a custom or more rigid AFO
- » Full range of sizes
- » Low-profile, effective for most indications
- » Inexpensive, labour savings for office or clinical application

Sizing Information

Size	Calf Circ (mm)	Foot Length (mm)	Brace Height (mm)	Code Left	Code Right
XSmall	152 - 190	101	165	B/3581	B/3586
Small	190 - 228	114	190	B/3582	B/3587
Medium	215 - 241	120	228	B/3583	B/3588
Large	228 - 226	133	254	B/3584	B/3589
XLarge	254 - 292	152	279	B/3585	B/3590