

ALL PURPOSE BOOT

The DARCO APB™ with rocker sole offers postoperative protection reducing plantar pressure under the forefoot and heel. The boot is commonly used in post-surgical healing, foot & ankle trauma and wounds or ulceration under metatarsal heads & toes.

A closed-toe design keeps dressings and toes warm and acts as a bumper when K-wires are present, providing additional room and comfort for the patient.



The perfect

SOLUTION

- » Used postoperatively for stability & mobilisation
- » Used with fibreglass cast bandages for the foot & ankle
- » Treatment of diabetic & insensitive foot



CONTACT US
TODAY!

E: info@beagleortho.com

T: 01254 268 788

W: beagleorthopaedic.com

Designed with great

FEATURES

- » Fully closed & water-repellent
- » Uses the same sole technology as the MedSurg healing shoes
- » Low-profile design: No height compensation necessary
- » Replaceable insoles available
- » The removable EVA insole is fully customisable to offload pressure points
- » The wider sole & deep opening provides unparalleled versatility by allowing for more room to accommodate large dressings, bandages & casts
- » Rocker bottom encourages a smooth roll through gait pattern & reduces plantar pressure on the heel & forefoot by over 25%
- » The wide ankle strap eliminates heel slippage as well as seats & secures the foot or cast in the shoe
- » The MetaShank™ construction provides a reinforced sole under the metatarsal heads to reduce motion
- » The strapless forefoot closure provides uniform security without the pressure of buckles or straps
- » Can be worn on the left or right foot

Sizing Information

Size	Shoe Size	Flat	Rocker Sole
XSmall	< 4	APQ0B	APB0B-ST
Small	4.5 - 6.5	APQ1B	APB1B-ST
Medium	7 - 8	APQ2B	APB2B-ST
Large	8.5 - 10	APQ3B	APB3B-ST
XLarge	10.5 - 12.5	APQ4B	APB4B-ST

Follow us on...

 BEAGLE ORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLETALK