

HEELWEDGE OFF-LOADING SHOE

The DARCO HeelWedge[®] is clinically proven to offload pressure by over 26% from the heel. Shifting weight to the mid and forefoot promotes faster healing after surgery, trauma or when wounds or ulcerations are present on the heel.



The Perfect

SOLUTION

- » Treatment of injuries or ulcers around the heel (heel ulcers)
- » Postoperative care after procedures involving plantar fasciitis



CONTACT US
TODAY!

E: info@beagleortho.com
T: 01254 268 788
W: beagleorthopaedic.com

Designed with great

FEATURES

- » Clinically proven to offload pressure by over 26% from the heel
- » Shifting weight to the mid and forefoot to promote faster healing after surgery, trauma or when wounds or ulcerations are present on the heel
- » Designed with a modified rocker bottom which encourages a smooth toe-off while maintaining a stable mid-foot base
- » Cloth toe protector keeps the toes warm and dressings clean & dry
- » Square toe design acts as a bumper to protect the toes providing a better fit for the left or right foot by eliminating toe overhang
- » Ankle strap seats the foot firmly into the shoe to reduce heel slippage & eliminate friction
- » The zoned outsole has an aggressive tread under the mid-foot where traction is needed most
- » The insole has twice the padding of standard insoles and can be modified
- » The insole can be completely removed to accommodate a DARCO Peg Assist customisable off-loading insole for an even more targeted off-load


Sizing Information

Size	Shoe Size	Code
XSmall	< 4	HQ0B
Small	4.5 - 6	HQ1B
Medium	6.5 - 8	HQ2B
Large	8.5 - 10	HQ3B
XLarge	10.5 - 12.5	HQ4B

Follow us on...

 BEAGLE ORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLETALK