

# ORTHOWEDGE OFF-LOADING SHOE

The DARCO OrthoWedge® is designed to protect the forefoot and promote faster healing after surgery, trauma, ulcers & infections. By reducing weight and pressure, a 57% reduction in metatarsal head pressure and 75% reduction in hallux toe pressure is achieved during the weight-bearing portion of the gait cycle.



*The perfect*

## SOLUTION

- » Used postoperatively after procedures to correct Hallux valgus, hammer toes/tailor's bunions
- » Forefoot ulcerations



CONTACT US  
**TODAY!**

**E:** [info@beagleortho.com](mailto:info@beagleortho.com)  
**T:** 01254 268 788  
**W:** [beagleorthopaedic.com](http://beagleorthopaedic.com)

*Designed with great*

## FEATURES

- » By reducing weight and pressure, a 57% reduction in metatarsal head pressure and 75% reduction in hallux toe pressure is achieved during the weight-bearing portion of the gait cycle
- » Designed with a 15° wedge sole, shifts weight back to the midfoot and heel
- » Square toe design acts as a bumper to protect the toes when K-wires are present. Providing a better fit for the left / right foot by eliminating toe overhang while the wide base provides extra stability
- » Ankle strap seats the foot firmly into the shoe reducing heel slippage and eliminating friction
- » The removable forefoot closure provides security eliminating buckle and strap pressure
- » The zoned outsole has an aggressive tread under the mid-foot where traction is needed most
- » The insole has twice the padding of standard insoles and can be modified
- » The insole can be completely removed to accommodate a DARCO Peg Assist customisable off-loading insole for an even more targeted off-load

### Sizing Information

Size	Shoe Size Europe	Shoe Size UK Men	Shoe Size UK Women	Code
XSmall	33.5 - 37	1.5 - 4.5	1 - 4	OQ0B
Small	37.5 - 39.5	4.5 - 6.5	4.5 - 6	OQ1B
Medium	40 - 41.5	7 - 7.5	6.5 - 7	OQ2B
Large	42 - 44	8 - 9.5	7.5 - 9.5	OQ3B
XLarge	44.5 - 47	10 - 12.5	10 - 12	OQ4B
Paediatric	Up to 31	12	12	PED-OW

## Follow us on...

 BEAGLE ORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLETALK