

Paediatric Shoulder Lok Brace

Material - Velfoam or Cotton



Manufactured
in the UK



Hand Wash &
Air Dry Only



ISO:13485
Accredited

Manufactured in the UK by Beagle manufacturing, the Shoulder Lok Brace is made from durable and lightweight material. The hook and loop fastening system is made from trimmable strapping and comprises moveable wide tips that allow for easy adjustments to suit the height and size of the Child.

The arm sling is designed to be comfortable to wear and easy to apply. The brace maintains the arm at 90 degrees and helps immobilise and support the shoulder, arm and wrist after injuries or surgery. The large pocket design allows room for a cast to aid in recovery and pain management. The sling aids recovery, by preventing key movements and stabilising the joints allowing the injured area to heal and providing pain relief, support and comfort. The sling can be used on either the left and right sides and can be made in a variety of colours and patterns. Contact customer services with a size and colour/pattern to order.



Indications for use

- » Immobilisation following injury or surgery

Features and benefits

- » Designed to immobilise & support the arm
- » Durable & lightweight material
- » Easy to apply
- » Comfortable to wear
- » Trimmable straps allow for easy adjustments
- » Moveable wide tips allow for a more precise fit
- » Universal: left or right
- » Can be hand washed & air dried



Sizing Information

| Size | Sling Width (mm) | Sling Height (mm) | Code | Colour / Pattern |
|--------|------------------|-------------------|----------|---|
| Small | 200mm | 150mm | SLP1/SLS | Specify from swatch & contact customer services |
| Medium | 230mm | 210mm | SLP2/SLS | Specify from swatch & contact customer services |
| Large | 280mm | 230mm | SLP3/SLS | Specify from swatch & contact customer services |

Any product supplied by Beagle Orthopaedic is intended for single patient use only. Any reuses of any product will invalidate the CE mark and any warranty.

SS/PD/023 REVO Date Issued: 04/10/2021

Beagle
ORTHOPAEDIC

Prospect Building, Duttons Way, Shadsworth Business Park, Blackburn, BB1 2QR
E: info@beagleortho.com W: beagleortho.com T: 01254 268 788