

Plantar Fasciitis Splint

Material - Black Elastic



Manufactured
in the UK



Wash Hand &
Air Dry Only



ISO:13485
Accredited

The Plantar Fasciitis splint is manufactured in the UK by Beagle Orthopaedic and is made from a durable black elastic. Comfortable to wear, the splint helps support the arch of your foot and provide gentle stretch at night or when resting, to reduce pain.

If tension and stress on the bowstring, in the arch of the foot, becomes too great, small tears can occur in the plantar fascia. These tears can become inflamed, causing pain commonly known as Plantar Fasciitis.



Indications for use

- » Plantar Fasciitis

Features and benefits

- » Durable material
- » Comfortable to wear, promoting a good stretch in the foot
- » The clever design uses an ankle strap, allowing for focused tension to be applied on the plantar fascia when the knee is in a bent or extended position
- » The superior strap (orange strips) allows the user to monitor pressure and progress
- » The semi-elastic inner sole gives flexibility to the plantar area and a firm hold at the toes, where it's needed. This allows the plantar fascia to stretch fully from the toes to the heel
- » The horizontal strap keeps your foot in place and contains a handy lower section in which you can place a Kool Pak to ease any inflammation or swelling that may have occurred

Sizing Information

Size	UK Shoe Size	Code
Small	3 - 7	PFS.SM
Medium	7.5 - 11	PFS.MD
Large	11.5 - 13	PFS.LG