

Orthomerica Supra-Lite™ AFO



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop. Injection moulded polypropylene allows for different thicknesses throughout the orthosis, thicker on the vertical aspect for rigidity, thinner footplate for ease of trimming, comfort and flexibility. It can be heat modified for minor fitting adjustments. The dorsum strap is included. It has removable / washable soft laminate calf padding for improved comfort and hygiene.



Indications for use

- » Flaccid foot drop

Features and benefits

- » Comfortable, low profile prefabricated leaf spring design to correct flaccid foot drop
- » Injection moulded polypropylene allows for different thickness throughout the orthosis, thicker on the vertical aspect for rigidity, thinner footplate for ease of trimming, comfort and flexibility
- » Can be heat modified for minor fitting adjustments
- » Dorsum strap included
- » Removable/washable soft laminate calf padding for improved comfort and hygiene
- » Can be hand washed and air dried only

Size	Male shoe size	Female shoe size	Footplate Length (mm)	Height (mm)	Left Code	Right Code
Small	< 6.5	< 6	241	266	B/3500	B/3505
Medium	7 - 9.5	6.5 >	280	292	B/3501	B/3506
Large	10 - 12.5	-	280	342	B/3502	B/3507
XLarge	13 >	-	298	393	B/3503	B/3508