

# BIOSKIN OA SPIRAL

## Knee Brace

Offloads the medial or lateral compartment to reduce osteoarthritis or meniscus pain

*The perfect*

### SOLUTION

- » Osteoarthritis
- » Meniscus Damage
- » Patellafemoral OA Pain

*Designed with great*

### FEATURES

- » The OA Spiral uses a three-point leverage system & a non-elastic strap to off-load the medial or lateral compartment to reduce pain
- » Magnetic buckle attachment on the spiral strap allows for quick & easy application & removal of the brace
- » Optional Q strap clinically proven to decrease pain in the PF joint by supporting patella alignment
- » Inner sleeve provides mild compression & holds the brace in place to allow easy application





# Clinically validated pain relief & high patient satisfaction!

## Lateral

### Circumference with knee fully extended

Size	Mid - patella	180mm above	Left Code	Right Code
Small	330 - 360	460 - 510	51501L-LAT	51501R-LAT
Medium	360 - 410	510 - 560	51502L-LAT	51502R-LAT
Large	410 - 430	560 - 640	51503L-LAT	51503R-LAT
X-Large	430 - 460	640 - 710	51504L-LAT	51504R-LAT
XX-Large	460 - 540	710 - 810	51505L-LAT	51505R-LAT

## Medial

### Circumference with knee fully extended

Size	Mid - patella	180mm above	Left Code	Right Code
Small	330 - 360	460 - 510	51501L-MED	51501R-MED
Medium	360 - 410	510 - 560	51502L-MED	51502R-MED
Large	410 - 430	560 - 640	51503L-MED	51503R-MED
X-Large	430 - 460	640 - 710	51504L-MED	51504R-MED
XX-Large	460 - 540	710 - 810	51505L-MED	51505R-MED

**LATEX FREE**  
**BIOMECHANICAL IMPROVEMENTS**  
**CONTOURED CUT**  
**ENHANCED PAIN RELIEF**  
**ADJUSTABLE SIZING**

CONTACT US **TODAY**

E: [info@beagleortho.com](mailto:info@beagleortho.com)

T: 01254 268 788

W: [beagleorthopaedic.com](http://beagleorthopaedic.com)

