

# DARCO

## ALLROUND SHOE

Beagle<sup>®</sup>  
ORTHOPAEDIC

The DARCO 'All-in-one' comprehensive healing, relief & rehabilitation shoe



### The perfect solution for:

- » Post-operative care
- » Wound care
- » Arthritis
- » Foot in risk category 0, 1 or 2

The AllRound healing shoe was designed to aid in the mobility of individuals with wide-ranging concerns & those with specific requirements, from Post-OP to wound care

CONTACT US  
TODAY!

E: [info@beagleortho.com](mailto:info@beagleortho.com)

T: 01254 268 788

W: [beagleorthopaedic.com](http://beagleorthopaedic.com)



MANUFACTURED  
IN LANCASHIRE 

**It's design supports comfort and ease-of-use to aid in the mobility of individuals with wide-ranging concerns and those with more specific requirements, with all the features that matter most, from Post-OP to wound care.**

*Designed with great*

## FEATURES:

- » The perfect, versatile option with a wide opening for easy access & plenty of room for bandages
- » The slip-resistant sole allows for easy mobility, offers good stability and shock absorption due to two different shore densities combining firmness and a cushioned sole for a softer landing
- » The shoe upper is designed in a light-weight and breathable nylon fabric that repels water & soil
- » Multi-layered closure for protection & a secure fit
- » A long ridged internal plate promotes weight distribution and the low profile design means there is no height compensation necessary
- » Can be worn on the right or left foot.

### Sizing Information

Size	Shoe Size	Code
XSmall	4 - 5.5	ARS0
Small	6 - 7	ARS1
Medium	7.5 - 8.5	ARS2
Large	9 - 10	ARS3
XLarge	10.5 - 11.5	ARS4
XXLarge	12 - 13	ARS5