

Healing Shoe

Features a rocker sole proven to reduce pressure on the forefoot and heel by 25%



1. BREATHABLE, LIGHTWEIGHT MESH

- ✓ Conforms to bony abnormalities of the foot
- ✓ Reinforced to stabilise the rear foot and reduce slippage

2. SQUARE TOE DESIGN

- ✓ Acts as a bumper to protect the toes when K-wires are present

3. ADJUSTABLE HOOK AND LOOP FASTENER

- ✓ Can fit up to a 17-inch circumference for large bandages and dressings
- ✓ The fastening system can be rotated to provide a left or right overlap and the ankle strap eliminates heel slippage

CONTACT US **TODAY**

E: info@beagleortho.com

T: 01254 268 788

W: beagleorthopaedic.com

The perfect

SOLUTION

- » Forefoot ulcerations
- » Post surgical healing and forefoot trauma
- » Postoperatively for stability and off-loading
- » Postoperatively after procedures to correct hallux valgus and hammer toes or tailor's bunions



**Patented MetaShank® forefoot support
which provides more rigid control under the
metatarsal heads and better protection**

Size	Male Shoe Size	Male Code	Female Shoe Size	Female Code
Small	< 8	MQM1B	< 4	MQW1B
Medium	8.5 - 10	MQM2B	4.5 - 5.5	MQW2B
Large	10.5 - 12	MQM3B	6 - 7.5	MQW3B
X-Large	12.5 - 14	MQM4B	-	-

MALE OR FEMALE SIZES
BREATHABLE
SOLID PROTECTION
COMFORTABLE
UNIVERSAL LEFT OR RIGHT

Follow us on...

 BEAGLE ORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLETALK