

# Shoulder Lok Sling

Beagle  
ORTHOPAEDIC

**An adjustable sling constructed in a soft, durable & lightweight Velfoam material that is bound at the edges**



The innovative strap design prevents pressure on the neck & supports the shoulder by maintaining position

**MANUFACTURED**  
IN LANGASHIRE 

## Designed with great Features:

- » Immobilises & supports the arm
- » Distributes weight evenly over the shoulders
- » Easy to apply & comfortable
- » Trimmable straps & movable wide tips allow for easy adjustments & a custom fit
- » Can be hand washed & air dried



**The adjustable straps distribute weight over the shoulders for enhanced comfort!**

## Indications for use:

- » Post-trauma
- » Post-operative
- » Immobilisation following injury or surgery

# CONTACT US TODAY!

**E:** [info@beagleortho.com](mailto:info@beagleortho.com)

**T:** 01254 268 788

**W:** [beagleorthopaedic.com](http://beagleorthopaedic.com)

**The sling aids recovery by preventing key movements & stabilising the joints allowing the injured area to heal & providing pain relief, support & comfort.**

### Sizing Information

Size	Sling Length (mm)	Sling Height (mm)	Code
Small	280mm	230mm	SL03/SLS
Medium	290mm	240mm	SL05/SLS
Large	330mm	255mm	SL07/SLS
XLarge	355mm	280mm	SL08/SLS

**FITS EASILY & SECURELY**  
**STABILISES THE ARM**  
**COMFORTABLE**  
**UNIVERSAL DESIGN**  
**COST-EFFECTIVE**



Also available  
in Paediatric

Contact our customer  
service team for more  
information!  
01254268788