

SPORLASTIC NEURO-LUX® II Shoulder Brace

A functional orthosis for stabilising the shoulder joint in cases of flaccid paryses after stroke

Designed with great Features:

- » Orthosis can be taken on and off with one hand
- » Different coloured plastic buckles make it easy to put on correctly
- » 2 lifting string: cross-joints, dynamic, for resetting force with the outward rotation on the shoulder joint
- » Breast belt: simple handling by closing and opening

**THERMOPLASTIC
MATERIAL PROMOTES
CORRECTION & AN
OPTIMAL FIT**

**CONTACT US
TODAY**

E: info@beagleortho.com

T: 01254 268 788

W: beagleorthopaedic.com

The perfect SOLUTION

- » Conservative therapy of the shoulder
- » Cranio-cerebral trauma
- » Hemiplegia
- » Paresis, flaccid paralysis of the shoulder and arm muscles
- » Peripheral nerve damage
- » Plexus injuries
- » Shoulder subluxation
- » Shoulder-hand syndrome
- » Trauma to the cervical spine



**FITS EASILY & SECURELY STABILISES THE ARM
COMFORTABLE SOFT & BREATHABLE MATERIAL**

Sizing Information

Size	Circ of lower arm (mm)	Chest Circ (mm)	Right	Left
1	210 - 250	750 - 900	07261/1R	07261/1L
2	210 - 250	900 - 1050	07261/2R	07261/2L
3	210 - 250	1050 - 1250	07261/3R	07261/3L
4	250 - 290	750 - 900	07261/4R	07261/4L
5	250 - 290	900 - 1050	07261/5R	07261/5L
6	250 - 290	1050 - 1250	07261/6R	07261/6L
7	290 - 350	750 - 900	07261/7R	07261/7L
8	290 - 350	900 - 1050	07261/8R	07261/8L
9	290 - 350	1050 - 1250	07261/9R	07261/9L