

# Standard Economy Wrist Brace 7" or 9"

Our Standard Economy Wrist Brace is made from a single weave elastic, providing support & comfort to the whole wrist region

## The Perfect Solution For...

- » Healing Fractures
- » Sprains
- » Repetitive strain injuries
- » Carpal tunnel syndrome

**Choose between 7"  
or 9" for varying  
levels of support for  
your wrist!**



## Great Features & Benefits:

- » Elastic provides support to the whole wrist region for added reassurance
- » Maintains full finger movement for rehabilitation exercises
- » Reduces pain by relieving strain on tendons
- » Can be hand wash & air dried

**A removable & easy to modify metal stay keeps the wrist in slight extension relieving strain on the tendons.**

## 7" WRIST BRACE

Sizing Information

Size	Wrist Circ (mm)	Code Box of 10		Code Singles	
		Left	Right	Left	Right
Small	130 - 150	EWB/013L	EWB/013R	EWB013L/SGL	EWB013R/SGL
Medium	150 - 170	EWB/015L	EWB/015R	EWB015L/SGL	EWB015R/SGL
Large	170 - 190	EWB/017L	EWB/017R	EWB017L/SGL	EWB017R/SGL
X-Large	190 - 210	EWB/018L	EWB/018R	EWB018L/SGL	EWB018R/SGL

## 9" WRIST BRACE

Sizing Information

Size	Wrist Circ (mm)	Code Box of 10		Code Singles	
		Left	Right	Left	Right
Small	130 - 150	EWB/913L	EWB/913R	EWB.SGL/93L	EWB.SGL/93R
Medium	150 - 170	EWB/915L	EWB/915R	EWB.SGL/95L	EWB.SGL/95R
Large	170 - 190	EWB/917L	EWB/917R	EWB.SGL/97L	EWB.SGL/97R
X-Large	190 - 210	EWB/918L	EWB/918R	EWB.SGL/98L	EWB.SGL/98R

**Maintain full finger movement for rehabilitation exercises**

**CONTACT US  
TODAY!**

**E:** [info@beagleortho.com](mailto:info@beagleortho.com)

**T:** 01254 268 788

**W:** [beagleorthopaedic.com](http://beagleorthopaedic.com)