

UNIVERSAL CLASBY

Humeral Brace

**Maintains skin integrity, aids fracture healing times
& improves clinical outcomes**



The perfect

SOLUTION

- » Mid & proximal diaphyseal fractures

Designed with great

FEATURES

- » Universal fit - wear it on the left or right side
- » The chest strap distributes weight across the body whilst not impinging on the breast area
- » Providing mild abduction to control & manage pain
- » There are 6 easy to remove metal stays that exert pressure equally to stabilise the fracture
- » 2 dual pull compression straps to apply compression and prevent rotation

Also available in:



Left and Right options and Paediatric

Beagle
ORTHOPAEDIC

CONTACT US **TODAY**

E: info@beagleortho.com

T: 01254 268 788

W: beagleorthopaedic.com

Clinically proven to control and manage pain!



Easy to apply, enhancing compliance and promoting early rehabilitation

Size	AC joint to elbow (mm)	Metal Stays Code	Plastic Stays Code
X-Small	330 - 360	UBC/12	UBC/22
Small	360 - 410	UBC/13	UBC/23
Medium	410 - 430	UBC/15	UBC/25
Large	430 - 460	UBC/17	UBC/27
X-Large	460 - 540	UBC/18	UBC/28

UNIVERSAL FIT
BREATHABLE
TRIMMABLE FOR A CUSTOM FIT
ANTIBACTERIAL
FULL CONTACT SUPPORT

Follow us on...

 BEAGLE ORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLEORTHOPAEDIC

 BEAGLETALK