

A functional ankle support that provides an alternative solution to tape bandaging



Designed with great FEATURES

- » Special knitted fabric on the instep for a wrinkle-free fit during rolling of the foot
- » Friction inserts made of silicone with raised nodules increase proprioception

Significant BENEFITS

- » Provides protection and easy to apply
- » 8 reins and 2 additional correction reins in the forefoot area, for improved foot elevation and pronation position through greater leverage



CONTACT US **TODAY!**

E: info@beagleortho.com
T: 01254 268 788
W: beagleorthopaedic.com

Figure 8 strap fastening stabilise the anterior talofibular ligament and lateral fastenings the calcaneofibular ligament

Indications for use:

- Primary and secondary prevention of sprains, especially in sports
- Treatment for acute stretched ligaments and chronic instability of the ankle

Sizing Information

Size	Above Ankle Circ (mm)	Left Code	Right Code
1	< 190	070731/L	070731/R
2	190 - 210	070732/L	070732/R
3	210 - 230	070733/L	070733/R
4	230 - 250	070734/L	070734/R
5	250 - 270	070735/L	070735/R

PROVIDES A SECURE FIT
EASY TO APPLY
COMFORTABLE TO WEAR
REINFORCES PROPRIOCEPTION
GENTLE ON SKIN

Follow us on...