

SPORLASTIC

Achillodyn Bandage

Provides proprioceptive muscular stabilisation and targeted compression



THE IDEAL SOLUTION FOR:

- » Achillodynia (tendinopathy) with/without partial rupture
- » Bursitis calcanea
- » Bursitis subachillea
- » Paratendinitis, Paratendinopathy
- » Postoperative foot and ankle

DESIGNED WITH GREAT FEATURES:

- » Proprioceptive muscular stabilisation
- » Separate heel wedge to temporarily reduce tension
- » Additional heel wedge for the healthy side to avoid scoliotic pelvis
- » The enhanced stretch makes the support easy to pull on, even when swelling is present



CONTACT US **TODAY**

Beagle®
ORTHOPAEDIC

E: info@beagleortho.com

T: 01254 268 788

W: beagleorthopaedic.com

ACHILLODYN

Septate pads gently massage the oedema away from the paratenon during movement and encourage blood flow

Sizing Information

Size	Above Ankle Circ (mm)	Code
1	< 190	0770711
2	190 - 210	0770712
3	210 - 230	0770713
4	230 - 250	0770714
5	250 - 270	0770715

The brace is constructed from an intelligent fabric that contains high absorption properties, meaning there is no floating effect experienced with full silicone supports

STRETCHY MATERIAL
INCREASED STABILITY
COMFORTABLE TO WEAR
MASSAGE EFFECT
HIGH COMFORT WEARING

Follow us on...