

SWEDE-O-UNIVERSAL

Ankle Brace

The ideal brace for the prevention and treatment of sprained ankles, providing optimal support during major physical stresses of the ankle joint



THE IDEAL SOLUTION FOR:

- » Heel pain
- » Plantar fasciitis
- » Ligament rupture and instability
- » Conservative therapy for sprained ankles
- » Postoperative follow-up treatment

Designed with lacing and lateral reinforcement, the offset panel traps the laces between the inner and outer layers to maintain the position for long periods

CONTACT US **TODAY**

Beagle®
ORTHOPAEDIC

E: info@beagleortho.com

T: 01254 268 788

W: beagleorthopaedic.com



DESIGNED WITH GREAT FEATURES:

- » Arch fits the contour of the foot
- » Internal U-shaped spiral stays provide extra support
- » Exclusive close spaced eyelets concentrate the holding power where it's needed most
- » Full elastic back ensures complete unrestricted blood flow

EFFECTIVE HEEL LOCK
SEAMLESS FIT
CONCENTRATED HOLDING POWER
UNRESTRICTED BLOOD FLOW
SUITABLE FOR USE WITH SHOES

Sizing Information

Size	Ankle Circ (mm)	Code
2	< 220	070722
3	220 - 240	070723
4	240 - 260	070724
5	260 - 280	070725
6	280 - 310	070726