

Orthomerica Carbon Supra-Lite™ AFO



Hand Wash &
Air Dry Only



Same Day
Dispatch



ISO:13485
Accredited

Low profile, dorsiflexion assist design. Increased rigidity in vertical aspect. This injection-moulded, carbon-reinforced polypropylene AFO fits the majority of patients requiring an AFO for Flaccid foot drop. Lightweight, comfortable design with carbon rebound feature provides dorsiflexion assist and the low arch and open heel allow it to fit in most shoes. Padded tibial crest strap and dorsum straps to enable accurate and secure positioning.



Indications for use

- » Flaccid foot drop
- » Post Stroke patients

Features and benefits

- » This new injection moulded, 15% carbon-reinforced polypropylene AFO has been designed to fit patients requiring an AFO for flaccid foot drop
- » Lightweight, comfortable design with carbon rebound feature provides dorsiflexion, and the open heels allow it to fit in most shoes
- » Padded tibial crest and dorsum straps enable accurate and secure positioning
- » Low profile, dorsiflexion assist design
- » Increased rigidity in vertical aspect
- » 4 sizes to fit the vast majority of patients

Sizing Information

Size	Footplate Length (mm)	Brace Height (mm)	Left	Right
Small	240	265	B/3520	B/3525
Medium	280	285	B/3521	B/3526
Large	280	330	B/3522	B/3527
XLarge	300	390	B/3523	B/3528